PAMAS
Partikelmess- und Analysesysteme



PAMAS S4031
Portable Particle Counting System for Liquids

Mobile and compact measuring instrument for raw water, waste water and potable water, for organic and corrosive liquids and for component cleanliness control

User-friendly operation using touch screen with graphic display

- The volumetric cell design of PAMAS sensors guarantees the highest accuracy, resolution and best statistical information
- Real portability with lab system accuracy
- Users can configure the system to their needs in profiles
- Highest repeatability and reproducibility
- Compact analysing system for liquids: raw water, waste water and potable water, organic and corrosive fluids
- Parts cleaning according to ISO 16232-9 and VDA-19
- Instrument can be used for batch and online sampling
- Sensor calibration with latex spheres according to ISO 21501
- Data storage of more than 4000 measurements

REV 09/2018

IN THE WORLD OF PARTICLES PAMAS COUNTS

PAMAS S4031

Portable particle counter with flexible software for various applications

PAMAS introduces a new development in general liquid monitoring. The **PAMAS S4031** particle counter is a portable self-contained instrument capable of analysing and reporting liquid contamination. The **PAMAS S4031** can be used for batch and online sampling.

Where instant results are essential the **PAMAS S4031** counter is the logical answer. Including our tried and tested high quality laser based sensor system, repeatable accurate results are produced quickly and easily.

With a simple-to-use touch screen user interface operation is easy and intuitive.

For higher concentration contamination measurements the instrument can be configured to measure up to 200,000 particles per ml. The unit has an integrated protection from contamination including a back flush operation to remove the contaminants from the system.

With 32 free adjustable size channels the system has the flexibility to match the most demanding applications.

The measuring results are reported according to common standards, including NAS 1638, SAE AS 4059 ISO 4406, GJB 420, GOST 17216 and NAVAIR 01-1A-17.


User-friendly download software for transfer of stored measurement data to a PC as a basic feature. Data files are compatible with most spreadsheet software.

TÜVRheinland ZERTIFIZIERT

Management System
ISO 9001:2015


www.tuv.com
ID 9105038017

Remote control with software PAMAS PMA



The population histogram shows cumulative and differential particle counts.

Automatic storage and documentation in readable format. Classification of particle number and size within the sensor's size range.



Complete documentation showing measuring results in absolute and relative numbers

Technical data

Sampling system:

- Wear resistant ceramic piston pump with controlled constant flow.

Particle Counter:

- 32-bit high performance CPU with sophisticated programmable digital domain signal conditioning and 4096 internal channels
- Data printout: 32 column thermo printer
- Data transfer: 8 bit ASCII code through USB port (57600 baud)
- Power supply: 90-230 V AC / 50-60 Hz or 12-30 V DC or via integrated battery for up to 3 hours operation
- Weight and Size: Approx. 8 kg 310 mm x 145 mm x 360 mm

Volumetric sensors:

PAMAS HCB-LD-50/50
Size range: 2 - 200µm
(1 - 400µm: option on request)
Max. particle concentration:
24,000 p/ml* at 25 ml/min**

PAMAS HCB-LD-25/25
Size range: 2 - 200 µm
Max. particle concentration:
120,000 p/ml* at 10 ml/min**

PAMAS HCB-LD-15/25
Size range: 2 - 100 µm
Max. particle concentration:
200,000 p/ml* at 10 ml/min**

Other particle sensors for larger particle sizes or higher particle concentrations are available on request.

* Coincidence error of 7.8%.
** Various flow-rates are available.

PAMAS HEAD OFFICE Dieselstraße 10, D-71277 Rutesheim, Phone: +49 7152 99 63 0, Fax: +49 7152 99 63-32, Email: info@pamas.de
PAMAS USA 1408 South Denver Avenue, Tulsa, OK 74119 USA, Phone: +1 918 743 6762, Fax: +1 918 743 6917, Email: clay.bielo@pamas.de
PAMAS BENELUX Mechelen Campus, Schaliënhoevedreef 20T, B-2800 Mechelen, Phone: +32 15 28 20 10, Mobile: +32 477 42 48 62, Email: paul.pollmann@pamas.de
PAMAS FRANCE Route du Tailleur 210/136, F-40170 Saint-Julien-en-Born, Mobile +33 6 25 33 20 41, Email: eric.colon@pamas.fr
PAMAS LATIN AMERICA Curitiba-Paraná, Brazil, Phone/Fax: +55 41 3022 5445, Mobile: +55 41 999 72 21 73, Email: marcelo.aiub@pamas.de
PAMAS INDIA No. 203, I floor, Oxford House, #15 Rustam Bagh Main Road, Bangalore 560017, India, Phone: +91 80 41 15 00 39, Email: info@pamas.in
PAMAS HISPANIA Calle Zubilleta No. 13 1°B, ES-48991 Algorta, Mobile: +34 67 75 39 699, Email: julian.malaina@pamas.de
PAMAS UK Sci-Tech Daresbury, Keckwick Lane, Daresbury, Cheshire WA4 4FS, Mobile: +44 79 17 71 33 66, Email: graeme.oakes@pamas.de

Please visit our website at www.pamas.de

® Registered trademarks are properties of their individual owners. All specifications are subject to change without notice.